

BC-PLAST Plastic Bottle Crates

Light is the *NEW STRONG* for Washing operations

1 Handy

Stainless Steel
Crate:
3,9 KG

BC-PLAST Plastic Crate:
2,4 KG

Mitigate the risk of burns post autoclavave with the plastic crates improved thermal insulation.



38%
lighter

Improved comfort

Less noise = less animal stress

Less friction = 38% reduction
in operator effort



3

Top Animal Welfare

Less friction, less noise when
removing plastic crate from the
trolley.

Simplified bottle release

Only 2,08N of force required to
slide the plastic crate from the
trolley (3,38N for the SS crate)

2

Ergonomic

Improved ergonomics:
access bottles without fully
removing the crate from
the trolley.

4

Long-lasting

The highly
durable plastic
crates are easy to
stack and do not
bend or distort
over time